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Development of Technique to Reduce Displacement of Focal Plane for X-ray Tube

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With prevalence of solid state detectors in CT scanner systems are becoming necessary various techniques to prevent degradation of image quality due to focal spot displacement of the x-ray tube unit. The already existing techniques are mechanisms to mechanically compensate focal spot displacement and correction algorithms applied on image data, but they have the disadvantages of cost-increase or longer reconstruction time. To solve this problem by the x-ray tube itself, the technique to minimize focal spot displacement has been developed.

Focal spot displacement is caused by swelling of the anode due to thermal expansion and flexion of the anode due to centrifugal force. In running x-ray tube, the rotor and the fixing part have higher temperature than the others. Then they have greater thermal expansion. We have conducted investigations in the fields of both structural mechanics and metallic material science and as the result, the high strength brazing technique and casting technique of compound structure have been established. In addition, the newly developed material was applied to the anode supporting component that expands to the opposite of the direction of the rotor expanded with heat. And it increased the correction effect.

The flexion of the anode due to high speed scanning is minimized with the new design including adoption of high-strength structure in a larger diameter of the anode and high-rigidity bearing for the rotor. The above development of the new technologies allows the focal spot displacement to be reduced 1/10 of the conventional displacement, thereby the compensation mechanism and correction algorithm are no more required and excellent image quality has been obtained by the CT system with the solid state detector.

Key Words: X-ray Tube, Displacement of Focal Plane, Ring Artifact

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