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Development of Acoustic Coupler for Elastography

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It is known that malignant diseases such as cancer generally become hardened, and therefore, Real-time Tissue Elastography* has allowed a judgment on whether a lesion is harder or softer than its surrounding tissues. The Strain Ratio, one of its measurement means, allows an objective judgment of the in vivo hardness by comparing the strain of a lesion with that of its surrounding tissues. This measurement requires a tissue as standard, but there may be a case where any appropriate standard tissue does not exist depending on the objective area. To solve this problem, we developed an acoustic coupler which can be a standard for the Strain Ratio. This paper reports on the Elastographic imaging method using an acoustic coupler together with some clinical cases.

Key Words: Elastography, Acoustic coupler, Strain Ratio

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- 1) Matsumura T, et al. : Development of Real-time Tissue Elastography. MEDIX 41, 30-35, 2004.
- 2) Ito A, et al. : Clinical application of US elastography for diagnosis. Radiology, 239(2) : 341-350, 2006.
- 3) Fukunari N : More Accurate and Sensitive Diagnosis of Thyroid Tumores with Elastography? Detection and Differential Diagnosis of Thyroid Cancers-. Medix Suppl. : 16-19, 2007.
- 4) Miyagawa T, et al. : Real-time Elastography for the Diagnosis of Prostate Cancer : Evaluation of Elastographic of Prostate Cancer : Acaluation of Elastographic Moving Images. Jpn J Clin Oncol, 39(6) : 394-398, 2009.
- 5) Tsutsumi M, et al. : Real-Time Balloon Inflation of Elastography for Prostate Cancer Detestion and Initial Evaluation of Clinicopatho ; pginalysis.AJR, 194 : W471-W476, 2010.
- 6) 藤本研治,ほか : Real-time Tissue Elastographyを用いた肝線維化の非侵襲的評価法～ Liver fibrosis index (LF Index)による stage判定. 肝臓, 51(9) : 539-541, 2010.
- 7) Tonomura A, et al. : Development of Strain Histogram Measurement Function and Clinical Applications in Hepatic Region. MEDIX 54, 37-41, 2011.
- 8) Waki K, et al. : Investigation of Strain Ratio Using Ultrasound Elastography Technique. Proc.ISICE : 449-452, 2007.
- 9) Ueno E, et al. : New quantitative method in brest elastography : Fat Lesion Ratio(FLR). Abstracts of RSNA 2007 ; LL-B2123-H04, 2007.
- 10) 植野映,ほか : 実践乳房超音波診断－基本操作,読影,最新テクニック,2 : 202-203 : 株式会社 山中書店, 2008.
- 11) Krouskop TA, et al. : Elastic Moduli of Breast and Prostate Tissue Under Compression. Ultrasonic Imaging, 20 : 260-274, 1998.