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Development of the FLR Assistance for the Strain Ratio Measurement in Breast Elastography

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We have developed a FLR Assistance. This is intended to improve the convenience and objectivity for one of the features of a Strain Ratio measurement of elastography. And that is what is measured semi-automatically strain ratio of a tumor area and a fat area in the mammary gland (FLR : Fat Lesion Ratio). In comparison with the value measured by the user manually, the average value is  $0.87 \pm 0.06SD$  in terms of the radius of the overlapping ratio of the tumor ROI and the fat ROI. We visually confirmed that the ROI on the fat layer of 93 cases in 99 cases. There was a strong positive correlation between the semi-automatic measurements and manual measurements ( $r = 0.84$ ,  $p < 0.0001$ ). Thus, this system is a promising tool to improve the throughput of Strain Ratio Measurement.

**Key Words:** FLR Assistance, Strain Ratio, Fat Lesion Ratio, FLR, Real-time Tissue Elastography, Ascendus

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- 1) Shiina T, et al. : Strain Imaging Using Combined RF and Envelope Autocorrelation Processing. Proc. of 1996 IEEE Ultrasonics Symp, 4 : 1331-1336, 1996.
- 2) Matsumura T, et al. : Development of Freehand Ultrasound Elasticity Imaging System and in vivo Results. First International Conference on the Ultrasonic Measurement and Imaging of Tissue Elasticity, 1 : 80, 2002.
- 3) Matsumura T, et al. : Development of Real-time Tissue Elastography. MEDIX 41, 30-35, 2004.
- 4) Matsumura T, et al. : Diagnostic results for breast disease by real-time elasticity imaging system. Proceedings of 2004 IEEE Ultrasonics Symposium : 1484-1487, 2004.

- 5) N.Nitta, T.Shiina : Estimation of Nonlinear Parameter of Tissues by Ultrasound. Japanese Journal of Applied Physics, 41, Part 1, 5B, 3572-3578, 2002.
- 6) Waki K, et al. : INVESTIGATION OF STRAIN RATIO USING ULTRASOUND ELASTOGRAPHY TECHNIQUE. Proc.ISICE 2007 p449-452
- 7) Itoh A, et al. : Clinical application of US elastography for diagnosis. Radiology, 239(2), 341-350, 2006.
- 8) Ueno E, et al. : New quantitative method in breast elastography : Fat Lesion Ratio(FLR). Abstracts of RSNA 2007 ; LL-BR2123-H04, 2007.