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Development of Shear Wave Measurement with a Reliability Index

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Ultrasound elastography starts from tumorous lesions in breast region using Strain elastography. In recent years, Shear wave elastography which detects the propagation of the shear wave velocity is widely used in chronic liver disease. We recently developed Shear Wave Measurement (SWM) which has a reliability index of the velocity measurement. Since SWM calculates the reliability index along with the shear wave velocity in every single measurement, it is possible to determine the appropriate measurement quantitatively. This function could contribute the improvement of the diagnostic performance. The purpose of this report is to describe a principle of SWM and reliability index, also show the elasticity phantom evaluation.

Key Words: Shear Wave Elastography, Reliability Index, Liver Phantom, Measurement

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

Methods Excitation methods	Strain imaging	Shear wave imaging
Manual compression	Strain elastography  RTE	N/A
Acoustic radiation force impulse excitation	ARFI Imaging 	Point shear wave speed measurement  SWM Shear wave speed imaging
Controlled external vibration		Transient elastography 

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- 1) Matsumura T, et al. : Development of Real-time Tissue Elastography. MEDIX, 41 : 30-35, 2004.
- 2) Shiina T : JSUM ultrasound elastography practice guidelines: basics and terminology. Journal of Medical Ultrasonics, 40 : 309-323, 2013.
- 3) Shiina T, et al. : WFUMB Guidelines and recommendations for clinical use of ultrasound elastography: Part 1: basic principles and terminology. Ultrasound Med. Biol., 41 : 1126-1147, 2015.
- 4) Nightingale K : Acoustic Radiation Force Impulse (ARFI) Imaging: A Review. Current Medical Imaging Reviews, 7 : 328-339, 2011.
- 5) Doherty J, et al. : Acoustic radiation force elasticity imaging in diagnostic ultrasound. IEEE Trans. Ultrason. Ferroelectr. Freq. Control., 60 : 685-701, 2013.

- 6) Kudo M, et al. : JSUM ultrasound elastography practice guidelines: liver. Journal of Medical Ultrasonics, 40 : 325-357, 2013.
- 7) Cohen-Bacrie C, et al. : QIBA Technical Committee for Shear Wave Speed (SWS) Measurement. RSNA 98th Scientific Assembly and Annu. Meeting (Poster), 2012.
- 8) Tonomura A, et al. : Development of strain histogram measurement function and clinical applications in hepatic region. MEDIX, 54 : 37-41, 2011.