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Development of a New Generation of Breast Real-time Tissue Elastography in Diagnostic Ultrasound System

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We have developed new generation Real-time Tissue Elastography\* (RTE) to improve objectivity and reproducibility on the breast elastography. The new generation RTE is consists of new scanning method, new Auto Frame Selection, new Assist Strain Ratio and zero-rejection on Strain Ratio, It can greatly remove the subjectivity of the operator by simple and easy operation.

On the phantom study, the measurement by the automatic combination resulted in low variability on inter/intra operator comparing to the measurement by manual.

It was shown that it has high reproducibility. New RTE is expected to be available with high reproducibility in clinical use.

**Key Words:** Strain Elastography, Auto Frame Selection, Assist Strain Ratio, No Manual Compression, Zero-rejection

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- 1) Cosgrove D, et al. : EFSUMB Guidelines and Recommendations on the Clinical Use of Ultrasound Elastography. Part 2: Clinical Applications.
- 2) Nakashima K, et al. : JSUM Ultrasound Elastography Practice Guidelines Breast. J Med Ultrasonics. 2013 ; 40 : 359-91.
- 3) Barr RG, et al. : WFUMB GUIDELINES AND RECOMMENDATIONS ON THE CLINICAL USE OF ULTRASOUND ELASTOGRAPHY : PART 2 ; BREAST
- 4) Nakashima K, et al. : Comprehensive ultrasound diagnosis for intraductal spread of primary breast cancer. Breast Cancer. 2013 ; 20 : 3-12.
- 5) Waki K, et al. : INVESTIGATION OF STRAIN RATIO USING ULTRASOUND ELASTOGRAPHY TECHNIQUE. In: Proceedings of ISICE 2007 : The First International Symposium on Information and Computer Elements 2007. 2007 ; 449-52.
- 6) Ueno E, et al. : New quantitative method in breast elastography: fat lesion ratio (FLR). [abstract]. In : Proceedings of the radiological society of North America scientific assembly and annual meeting. Oak Brook, IL: Radiological Society of North America. 2007 ; 697.
- 7) Umemoto T, et al. : Ex vivo and in vivo assessment of the non-linearity of elasticity properties of breast tissues for quantitative strain elastography. Ultrasound Med Biol. 2014 Aug ; 40(8) : 1755-68.
- 8) Baba H, et al. : Development of Assist Strain Ratio for the Strain Ratio Measurement in Breast Elastography. MEDIX. 2013 ; 58 : 42-5.
- 9) Ueno E, et al. : A Preliminary Prospective Study to Compare the Diagnostic Performance of Assist Strain Ratio versus Manual Strain Ratio. J Med Ultrasonics. DOI : 10.1007/s10396-015-0633-5.