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Function of Real-time Virtual Sonography in the period of multipolar radiofrequency ablation for patients with hepatocellular carcinoma

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Multipolar radiofrequency ablation (RFA) is expected to produce a large ablative zone and reduce intrahepatic recurrence of hepatocellular carcinoma (HCC) by no touch ablation. Nevertheless, it requires considerable skill and it is difficult to get the intended ablative zone. ARIETTA^{*2} 850, a new ultrasound system, bundles 3D Sim-Navigator and E-field simulator through advances in function of Real-time Virtual Sonography^{*1}. These new function could be contributed to decrease the difficulty in multipolar RFA.

Key Word: Multipolar Radiofrequency Ablation, Hepatocellular Carcinoma, 3D Sim-Navigator, E-field Simulator

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- 1) Obara K, et al. Insufficient radiofrequency ablation therapy may induce further malignant transformation of hepatocellular carcinoma. Hepatol Int, 2 :116-123, 2008.
- 2) Nakazawa T, et al. Radiofrequency ablation of hepatocellular carcinoma: correlation between local tumor progression after ablation and ablative margin. AJR Am J Roentgenol, 188:480-488, 2007
- 3) 大久保 裕直 肝がん:診断と治療の進歩 ラジオ波焼灼術. 日内誌 103:78-86, 2014
- 4) Shiina S, et al. Radiofrequency ablation for hepatocellular carcinoma: 10-years outcome and prognostic factors. Am J Gastroenterol,107:569-577, 2012
- 5) Koda M, et al. Complications of radiofrequency ablation for hepatocellular carcinoma in a multicenter study: An analysis of 16346 treated nodules in 13283 patients. Hepatol Res, 42 :1058-1064, 2012.
- 6) Imamura J, et al. Neoplastic seeding after radiofrequency ablation for hepatocellular carcinoma. Am J Gastroenterol,103 :3057-3062, 2008.
- 7) Masuda T, et al: Intrahepatic dissemination of hepatocellular carcinoma after local ablation therapy. J Hepatobiliary Pacreat Surg,15 :589-595, 2008.
- 8) Osaki Y, et al. Clinical effectiveness of bipolar radiofrequency ablation for small liver cancers. J Gastroenterol, 48: 874-883, 2013.
- 9) Lee JM, et al. A comparative experimental study of in-vivo efficiency of hypertonic saline-enhanced hepatic bipolar and monopolar radiofrequency ablation. Korean J Radiol, 4: 163-169, 2003.
- 10) Kawamura Y, et al. Potential of a no-touch pincer ablation procedure that uses a multipolar radiofrequency ablation system to prevent intrasubsegmental recurrence of small and single hepatocellular carcinomas. Hepatol Res,47:1008-1020, 2017
- 11) Hirooka M, et al. Usefulness of a New Three-Dimensional Simulator System for Radiofrequency Ablation. PLoS One, 11:e0148298, 2016