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Image Diagnosis Support using Hybrid Learning

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The use of AI(Artificial Intelligence) has been activated to solve problems such as improving work efficiency at medical sites, improving the quality of diagnosis, and improving medical outcome. Currently, most of the related study is still research level, but in the United States some medical devices that are approved by the FDA(Food and Drug Administration) are starting to appear, and movements toward practical use are also proceeding. Under the concept of “Hybrid Learning”, Hitachi is promoting R & D on original diagnostic imaging AI technology that leads to highly accurate and medical validity results. In this paper, as a specific example, we describe lung cancer CAD technology that realizes highly accurate lesion detection function and MRI high image quality technology to realize high speed imaging.

Key Word: AI, Deep Learning, Radiology, Diagnostic Support, Machine Learning

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
(Multi-Adaptive Convolutional Neural Network Reconstruction : MA-CNNR)

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Method	PSNR [dB]
3層CNN	32.8
3層CNN+Residual	33.8
MA-CNNR	35.1

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- 1) https://www.accessdata.fda.gov/cdrh_docs/pdf17/K170195.pdf
- 2) <https://www.fda.gov/newsevents/newsroom/pressannouncements/ucm604357.htm>
- 3) The National Lung Screening Trial Research Team : "Reduced lung-cancer mortality with low-dose computed tomographic screening", N. Eng. J. Med., 365, 395-409, 2011.
- 4) Kusano S, et al. : "Efficacy of computer-aided diagnosis in lung cancer screening with low-dose spiral computed tomography: receiver operating characteristic analysis of radiologists' performance", Jpn. J. Radiol, 28 (9), 649-655, 2010.
- 5) 影山昌広, ほか: "新しい畳み込みニューラルネットワーク設計法を用いた肺がんCT検診向け結節検出CADシステム", 第75回日本放射線技術学会総会学術大会, 2019
- 6) Armato III, Samuel G., McLennan, Geoffrey, Bidaut, Luc, McNitt-Gray, Michael F., Meyer, Charles R., Reeves, Anthony P., ... Clarke, Laurence P. (2015). Data From LIDC-IDRI. The Cancer Imaging Archive. <http://doi.org/10.7937/K9/TCIA.2015.LO9QL9SX>
- 7) D. L. Donoho : "Compressed sensing", IEEE Trans. Inf. Theory, vol.52, no.4, pp.1289-1306, April. 2006.
- 8) Chao Dong, et al. : "Learning a Deep Convolutional Network for Image Super-Resolution", ECCV, 2014.
- 9) Li Xu, et al. : "Deep Convolutional Neural Network for Image Deconvolution", NIPS, 2014
- 10) Jiwon Kim, et al. : "Accurate Image Super-Resolution Using Very Deep Convolutional Networks", arXiv: 1511.04587v1, 2015
- 11) 野口喜実, ほか : "クラスタリングを用いたCNN選択によるMRI高画質化の検討", JAMIT2018
- 12) IXI Dataset , brain-development.org, <http://brain-development.org/ixi-dataset/>